



**The Leader
in First Aid
Emergency
Treatment
for burns.**

SOFT-SIDED BURN KITS



HARD-SIDED BURN KITS



BE PREPARED FOR ANY BURN EMERGENCY.

More than 25 years ago, WATER-JEL introduced water-based gel technology that revolutionized the emergency first aid treatment of burns. Today, WATER-JEL is used by EMS and industrial health and safety professionals around the world.

The “immediate cooling” of burns is recommended by most EMS protocols and curriculum.

The gel used in WATER-JEL Fire Blankets and Dressings is made from 96% water and is very slow to evaporate, so blankets and dressings remain moist and won't stick to burn injuries. They provide thermal (not evaporative) cooling of burns without the need for re-wetting so they're easy to use and won't get too cold or increase the risk of hypothermia – even on large extent burns. The gel also contains medically pure tea tree oil, a natural antibacterial.

- Leading the World in Emergency Burn Care!
- Exclusive Supplier to the US Military
- Mandated in Formula One Racing . . . they cannot start a race without it trackside in every pit crew
- Made in the USA at an FDA-regulated facility
- Used by many Fortune 500 companies worldwide
- Helps reduce the impact of a burn and lost work time

WATER-JEL Burn Kits come in different product assortments for the different work environments of First Responders, Safety Managers and Medical Professionals.



it
COOL

WATER-JEL BURN KITS MEET THE NEEDS OF EVERY INDUSTRY.



WATER-JEL SOFT-SIDED EMERGENCY BURN KIT

#EBK2-3 - Large w/Fire Blanket

One 12" x 16" Burn Dressing (with pre-cut holes for facial burns)
 One 8" x 18" Burn Dressing
 One 4" x 16" Burn Dressing
 Three 4" x 4" Burn Dressing
 Four sterile gauze roll 3" x 4 yds.

One 36" x 30" (3' x 2.5') Fire Blanket/Burn Wrap
 One scissors
 One instructions card
 Heavy-duty nylon carry bag with inner pockets.

- Vital in areas where high heat flames, arc flash, hot liquid, steam or chemicals exist.
- Everything needed for Industrial, Institutional and Medical burn emergencies.

Case Pack

3 Kits
Case Size:
 15.38" x 12.63" x 8.13"
Case Weight:
 22.3 lbs
Case Cube:
 0.92 cu. ft.



WATER-JEL SOFT-SIDED EMERGENCY BURN KIT

#EBK1-6 - Small

One 12" x 16" Burn Dressing (with pre-cut holes for facial burns)
 Two 4" x 4" Burn Dressing
 Two 4" x 16" Burn Dressing
 Two sterile gauze roll 3" x 4 yds.

One scissors
 One instructions card
 Heavy-duty nylon carry bag with inner pockets.

- Perfect for use in work vehicles.
- Durable nylon handle straps for easy carrying, or can be hung on a wall.

Case Pack

6 Kits
Case Size:
 15.38" x 12.63" x 8.13"
Case Weight:
 13.7 lbs
Case Cube:
 0.92 cu. ft.



WATER-JEL INDUSTRIAL/WELDING BURN KIT

#IWK-5

One 12" x 16" Burn Dressing (with pre-cut holes for facial burns)
 One 8" x 18" Burn Dressing
 One 4" x 16" Burn Dressing
 Two 4" x 4" Burn Dressing

Two 2" x 6" Burn Dressing
 One BurnJel dispenser with 25 unit dose packets
 One instructions card • One elastic bandage

- Carry Case: High-impact polystyrene kit box has living hinges, strong nylon catches, moisture-resistant cover gasket and fold-down combination handle/hanger. Can be wall-mounted in locations where burns can happen.

Case Pack

5 Kits
Case Size:
 16" x 12" x 8"
Case Weight:
 23.3 lbs
Case Cube:
 0.89 cu. ft.



WATER-JEL FOOD SERVICES BURN KIT

#FSK-5

One 12" x 16" Burn Dressing (with pre-cut holes for facial burns)
 Two 4" x 4" Burn Dressing
 Two 2" x 6" Burn Dressing

One BurnJel dispenser with 25 unit dose packet
 One instructions card • One elastic bandage

- Carry Case: High-impact polystyrene kit box has living hinges, strong nylon catches, moisture-resistant cover gasket and fold-down combination handle/hanger. Can be wall-mounted in locations where burns can happen.

Case Pack

5 Kits
Case Size:
 16" x 12" x 8"
Case Weight:
 12.3 lbs
Case Cube:
 0.89 cu. ft.



WATER-JEL LARGE EMERGENCY BURN STATION

#EBSL-5

One Face Dressing 12" x 16"
 One Dressing 8" x 18"
 Two Dressings 4" x 16"
 Two Dressings 4" x 4"

One Burn Jel 2 oz bottle
 10 Burn Jel unit dose 3.5g
 One pair of scissors
 Two rolls of gauze

- Provides immediate access to a wide range of burn treatment products.
- Place in strategic locations where burn accidents are most likely to occur.

Case Pack

5 stations
Case Size:
 17.63" x 22" x 12.75"
Case Weight:
 12.55 lbs
Single Pack Size:
 17.25" x 21.25" x 2.25"
Weight:
 2.33 lbs



WATER-JEL SMALL EMERGENCY BURN STATION

#EBSS-5

Two Dressings 2" x 6"
 Two Dressings 4" x 4"

6 Burn Jel unit dose 3.5g
 One roll of gauze

- Offers full product visibility.
- Very well suited for fast food operations and smaller businesses.

Case Pack

5 stations
Case Size:
 15.5" x 8.75" x 9.25"
Case Weight:
 2.73 lbs
Single Pack Size:
 14.875" x 8.5" x 1.5"
Weight:
 .47 lbs



50 Broad Street, Carlstadt, NJ 07072
 Tel: (201) 507-8300 • Fax (201) 507-8325 • (800) 275-3433
 E-mail: info@waterjel.com • Website: www.waterjel.com